




Factory Address: 4# Workshop, No. 2875 Dongfu Avenue,  
Haicang District, Xiamen  
Tel.: +86 592 2293630  
Fax: +86 592 2296092  
eMail: sales@bellamoon.cn  
Website: <https://bellamoon.en.alibaba.com>

**倍乐梦（厦门）医疗科技有限公司**  
Bellamoon (Xiamen) Medical Technology Co.,Ltd



▶ Professional Manufacturer of Ward Nursing Equipment

倍乐梦（厦门）医疗科技有限公司  
Bellamoon (Xiamen) Medical Technology Co.,Ltd

## BC204R 2 Crank manual hospital bed

- Material: Cold rolled steel tube bed frame, ABS plastic bed board.
- Bed size: 215\*98cm (Length \* width)
- Board size: 196\*90cm (Length \* width), Height of bed above the ground: 52cm
- Basic parameters: Folding angle of back panel 0°~70° Folding angle of leg plate 0°~30°
- Configuration: Central control silent casters, detachable PP headboard and foot board, Aluminum alloy folding guardrail, infusion stand.



## BC205R medical clinic patient beds

- Material: Cold rolled steel pipe bed frame, steel plate bed surface, electrostatic powder coating on the surface.
- Bed size: length: 215cm, width: Guardrail raised: 105cm, Guardrail down: 99cm
- Board size: 196\*90cm (Length \* width), Height of bed surface: 50cm
- Basic parameters: Folding angle of back panel 0°~70°  
Folding angle of leg plate 0°~30°
- Configuration: Central control silent casters, detachable PP headboard, footboard, PP guardrail, Backplane angle indicator, infusion stand, anti-collision pad, 2 cranks.



## BC306R Adjustable 3 crank manual hospital bed

- Material: Cold rolled steel pipe bed frame, steel plate bed surface
- Bed size: 215\*98cm (length \* width)
- Board size: 196\*90cm (Length \* width), Height of bed above the ground: Highest: 70cm, Lowest: 45cm
- Basic parameters: Folding angle of back panel  $0^{\circ}\sim\geq 70^{\circ}$   
Folding angle of leg plate  $0^{\circ}\sim\geq 30^{\circ}$
- Configuration: Central control silent casters, detachable PP headboard, PP guardrail, infusion stand (optional)




## BC307R Manual multifunctional patient bed

- Material: Cold rolled steel pipe bed frame, steel plate bed surface, surface electrostatic powder spraying.
- Bed size: 215\*98cm (length \* width)
- Board size: 196\*90cm (Length \* width), Height of bed surface: Highest: 70cm  
Lowest: 45cm
- Basic parameters: Folding angle of back panel  $0^{\circ}\geq 70^{\circ}$   
Folding angle of leg plate  $0^{\circ}\geq 30^{\circ}$
- Configuration: Silent casters, detachable PP headboard, aluminum alloy folding guardrail, infusion stand (optional)



## BC2004A 2 Functions Manual medical bed

- Material: Cold rolled steel tube bed frame.
- Bed size:2060\*970\*530(Length\*width\*height)
- Board size:196\*90cm(Length \* width),Height of bed above the ground: 52cm
- Basic parameters:Folding angle of back panel 0°~≥70°  
Folding angle of leg plate 0°~≥30°
- Configuration:L:4 independent brake silent casters, detachable PP headboard and foot board,Aluminum alloy folding guardrail, 1 platform shoe rack, infusion stand.



## BC2005A hospital bed adjustable patient bed

- Material: Cold rolled steel pipe bed frame, steel plate bed surface, electrostatic powder coating on the surface.
- Bed size:length:208\*98\*47.5cm(Length\*width\*height)
- Board size:196\*90cm(Length \* width),Height of bed surface:48cm
- Basic parameters:Folding angle of back panel 0°~≥70°  
Folding angle of leg plate 0°~≥30°
- Configuration:4 independent brake silent casters, detachable PP headboard, foot board, aluminum alloy folding guardrail, 1 platform shoe rack, 1 infusion stand, 2 double cranks, 4 anti-collision cushions.



## BC202R 2 Crank luxury hospital bed

- Material: Cold rolled steel pipe bed frame, steel plate bed surface
- Bed size: Length: 215cm, Width: when the guardrail is raised: 105cm, down: 99cm.
- Board size: 196\*90cm (Length \* width), Height of bed above the ground: 52cm.
- Basic parameters: Folding angle of back panel  $0^{\circ} \sim \geq 70^{\circ}$   
Folding angle of leg plate  $0^{\circ} \sim \geq 30^{\circ}$
- Configuration: Central control silent casters, detachable PP headboard, PP guardrail, Backplane angle indicator, infusion stand (optional)



## BC401R 4 Crank medical luxury hospital bed

- Material: Cold rolled steel pipe bed frame, steel plate bed surface, surface electrostatic powder spraying.
- Bed size: 223cm, Width: when the guardrail is raised: 105cm, down: 99cm.
- Board size: 196\*90cm (Length \* width), Height of bed surface: Highest: 70cm  
Lowest: 45cm
- Basic parameters: Folding angle of back panel  $0^{\circ} \geq 70^{\circ}$   
Folding angle of leg plate  $0^{\circ} \geq 30^{\circ}$   
The overall forward inclination angle of the bed surface  $0^{\circ} \sim \geq 12^{\circ}$   
The whole bed surface tilt angle  $0^{\circ} \sim \geq 12^{\circ}$
- Configuration: Silent casters, detachable PP headboard, aluminum alloy folding guardrail, infusion stand (optional)



## BC308R Three function manual medical bed

- Material: Cold rolled steel pipe bed frame, steel plate bed surface.
- Bed size: 215\*98cm (length \* width)
- Board size: 196\*90cm (Length \* width), Height of bed above the ground: Highest: 70cm, Lowest: 45cm
- Basic parameters: Folding angle of back panel  $0^{\circ} \sim \geq 70^{\circ}$   
Folding angle of leg plate  $0^{\circ} \sim \geq 30^{\circ}$
- Configuration: Silent casters, detachable PP headboard, aluminum alloy folding guardrail, infusion stand (optional)

